


# GeneSTAR

## StarPAGE Precast Gel



**GeneSTAR 공식대리점**

TEL[1] 02-981-4158

TEL[2] 010-5266-6018

E-mail. [withsbio@naver.com](mailto:withsbio@naver.com)



**위드에스바이오** 

카카오톡 채널 추가로 상담하기!

# Selection Guide

4-12% MOPS	4-20% MOPS	8% MOPS	10% MOPS	12% MOPS
270 Kd	270 Kd	270 Kd	185 Kd	185 Kd
185 Kd	185 Kd	185 Kd	140 Kd	140 Kd
140 Kd	140 Kd	140 Kd	115 Kd	115 Kd
115 Kd	115 Kd	115 Kd	80 Kd	80 Kd
80 Kd	80 Kd	115 Kd	65 Kd	65 Kd
65 Kd	65 Kd	80 Kd	50 Kd	50 Kd
50 Kd	50 Kd	65 Kd	40 Kd	40 Kd
40 Kd	40 Kd	50 Kd	30 Kd	30 Kd
30 Kd	30 Kd	40 Kd	25 Kd	25 Kd
25 Kd	15 Kd	30 Kd	15 Kd	15 Kd
15 Kd	10 Kd	25 Kd		10 Kd

## Ordering Information

Cat.No.	Product	Size	Price
GPG0811	StarPAGE™ 8% 11Wells	10/Box	90,000
GPG0815	StarPAGE™ 8% 15Wells	10/Box	90,000
GPG1011	StarPAGE™ 10% 11Wells	10/Box	90,000
GPG1015	StarPAGE™ 10% 15Wells	10/Box	90,000
GPG1211	StarPAGE™ 12% 11Wells	10/Box	90,000
GPG1215	StarPAGE™ 12% 15Wells	10/Box	90,000
GPG4111	StarPAGE™ 4~12% 11Wells	10/Box	90,000
GPG4115	StarPAGE™ 4~12% 15Wells	10/Box	90,000
GPG4211	StarPAGE™ 4~20% 11Wells	10/Box	90,000
GPG4215	StarPAGE™ 4~20% 15Wells	10/Box	90,000
GBM0080	MOPS-SDS Running Buffer (10X)	500mL	50,000
GFT0002	StarFast Transfer Buffer (20X)	500mL	150,000

(\* Note : Tris-Glycine electrophoresis buffer is not compatible with the Bis-Tris buffer system of StarPAGE Precast Gel, so do not use it.)



**GeneSTAR 공식대리점**  
 TEL[1] 02-981-4158  
 TEL[2] 010-5266-6018  
 E-mail. withsbio@naver.com



카카오톡 채널 추가로 상담하기!